

For Research Use Only

ACF1 Monoclonal antibody, PBS Only

Catalog Number: 67947-1-PBS



Basic Information

Catalog Number: 67947-1-PBS	GenBank Accession Number: NM_013448	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 11177	CloneNo.: 1C4F8
Source: Mouse	UNIPROT ID: Q9NRL2	
Isotype: IgG2b	Full Name: bromodomain adjacent to zinc finger domain, 1A	
Immunogen Catalog Number: AG30607	Calculated MW: 179 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67947-1-Ig (ACF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67947-1-PBS in a different storage buffer formulation.