

For Research Use Only

Endophilin 1 Monoclonal antibody, PBS Only



Catalog Number: 68029-1-PBS

Basic Information

Catalog Number: 68029-1-PBS	GenBank Accession Number: BC032825	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6456	CloneNo.: 1G4H1
Source: Mouse	UNIPROT ID: Q99962	
Isotype: IgG2a	Full Name: SH3-domain GRB2-like 2	
Immunogen Catalog Number: AG24793	Calculated MW: 39.8 kDa	
	Observed MW: 38-40 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat, pig, rabbit

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68029-1-Ig (SH3GL2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68029-1-PBS in a different storage buffer formulation.