

For Research Use Only

RABGAP1L Monoclonal antibody, PBS Only



Catalog Number: 68356-1-PBS

Basic Information

Catalog Number:

68356-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG5305

GenBank Accession Number:

BC041888

GeneID (NCBI):

9910

UNIPROT ID:

B7ZAP0

Full Name:

RAB GTPase activating protein 1-like

Calculated MW:

815 aa, 93 kDa

Observed MW:

100 kDa, 140 kDa

Purification Method:

Protein A purification

CloneNo.:

4A3A9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

RABGAP1L, also known as TBC1D18, is a member of the Tre2/Bub2/Cdc16 (TBC)-domain-containing protein family. Overexpression of RABGAP1L potentiated the antiviral action of IFN against a panel of viruses, including IAV, and its function required a catalytically active TBC domain.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68356-1-Ig (RABGAP1L antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68356-1-PBS in a different storage buffer formulation.