

For Research Use Only

RFC4 Monoclonal antibody, PBS Only

Catalog Number: 68390-1-PBS



Basic Information

Catalog Number: 68390-1-PBS	GenBank Accession Number: BC017452	Purification Method: Protein G purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 5984	CloneNo.: 1A6A11
Source: Mouse	UNIPROT ID: P35249	
Isotype: IgG1	Full Name: replication factor C (activator 1) 4, 37kDa	
Immunogen Catalog Number: AG33295	Calculated MW: 37 kDa	
	Observed MW: 37-40 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68390-1-Ig (RFC4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 68390-1-PBS in a different storage buffer formulation.