

For Research Use Only

HSDL2 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68414-3-PBS



Basic Information

Catalog Number: 68414-3-PBS	GenBank Accession Number: BC004331	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 84263	CloneNo.: 2H2E7
Source: Mouse	UNIPROT ID: Q6YN16	
Isotype: IgG1	Full Name: hydroxysteroid dehydrogenase like 2	
Immunogen Catalog Number: AG8314	Calculated MW: 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

68414-3-PBS targets HSDL2 as part of a matched antibody pair.

MP50441-1: 68414-2-PBS capture and 68414-3-PBS detection (validated in Cytometric bead array)

MP50441-2: 68414-3-PBS capture and 68414-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50441-1, HSDL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68414-2-PBS. Detection antibody: 68414-3-PBS. Standard: Ag8314. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50441-2, HSDL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68414-3-PBS. Detection antibody: 68414-4-PBS. Standard: Ag8314. Range: 0.098-100 ng/mL.