

For Research Use Only

STX4 Monoclonal antibody, PBS Only

Catalog Number: 68429-1-PBS



Basic Information

Catalog Number:

68429-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7053

GenBank Accession Number:

BC002436

GeneID (NCBI):

6810

UNIPROT ID:

Q12846

Full Name:

syntaxin 4

Calculated MW:

34 kDa

Observed MW:

35 kDa

Purification Method:

Protein G purification

CloneNo.:

4D6C3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68429-1-Ig (STX4 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68429-1-PBS in a different storage buffer formulation.