

For Research Use Only

Synaptotagmin-9 Monoclonal antibody, PBS Only



Catalog Number: 68430-1-PBS

Basic Information

Catalog Number: 68430-1-PBS	GenBank Accession Number: BC029605	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 143425	CloneNo.: 3D2F2
Source: Mouse	UNIPROT ID: Q86556	
Isotype: IgG1	Full Name: synaptotagmin IX	
Immunogen Catalog Number: AG33194	Calculated MW: 491 aa, 52 kDa	
	Observed MW: 52 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat, rabbit

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68430-1-Ig (Synaptotagmin-9 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68430-1-PBS in a different storage buffer formulation.