

For Research Use Only

USP9X Monoclonal antibody, PBS Only

Catalog Number: 68460-1-PBS

Featured Product



Basic Information

Catalog Number:

68460-1-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_001039590

GeneID (NCBI):

8239

UNIPROT ID:

Q93008

Full Name:

ubiquitin specific peptidase 9, X-linked

Calculated MW:

292 kDa

Observed MW:

260 kDa

Purification Method:

Protein G purification

CloneNo.:

1H6G11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68460-1-Ig (USP9X antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68460-1-PBS in a different storage buffer formulation.



WB result of USP9X antibody (68460-1-Ig; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-USP9X transfected HEK-293 cells. This data was developed using the same antibody clone with 68460-1-PBS in a different storage buffer formulation.