

For Research Use Only

UBXN1 Monoclonal antibody, PBS Only



Catalog Number: 68499-1-PBS

Basic Information

Catalog Number:

68499-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33606

GenBank Accession Number:

BC000902

GeneID (NCBI):

51035

UNIPROT ID:

Q04323

Full Name:

UBX domain protein 1

Calculated MW:

33 kDa

Observed MW:

40 kDa

Purification Method:

Protein G purification

CloneNo.:

3A6C2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68499-1-Ig (UBXL1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 68499-1-PBS in a different storage buffer formulation.