

For Research Use Only

MTHFD2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68524-2-PBS

Basic Information

Catalog Number: 68524-2-PBS	GenBank Accession Number: BC017054	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10797	CloneNo.: 1B7H2
Source: Mouse	UNIPROT ID: P13995	
Isotype: IgG1	Full Name: methylenetetrahydrofolate dehydrogenase (NADP+ dependent) 2, methenyltetrahydrofolate cyclohydrolase	
Immunogen Catalog Number: AG33605	Calculated MW: 350 aa, 38 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

68524-2-PBS targets MTHFD2 as part of a matched antibody pair:

MP50286-1: 68524-2-PBS capture and 68524-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50286-1, MTHFD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68524-2-PBS. Detection antibody: 68524-3-PBS. Standard: Ag33605. Range: 0.391-100 ng/mL.