

For Research Use Only

# UBA5 Monoclonal antibody, PBS Only

Catalog Number: 68537-1-PBS



## Basic Information

Catalog Number:

68537-1-PBS

Size:

100ug, Concentration: 1 mg/ml by  
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33685

GenBank Accession Number:

BC009737

GeneID (NCBI):

79876

UNIPROT ID:

Q9GZZ9

Full Name:

ubiquitin-like modifier activating  
enzyme 5

Calculated MW:

404 aa, 45 kDa

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

3E5C8

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68537-1-Ig (UBA5 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68537-1-PBS in a different storage buffer formulation.