

For Research Use Only

MSX2 Monoclonal antibody, PBS Only



Catalog Number: 68550-1-PBS

Basic Information

Catalog Number:

68550-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33220

GenBank Accession Number:

BC015509

GeneID (NCBI):

4488

UNIPROT ID:

P35548

Full Name:

msh homeobox 2

Calculated MW:

29 kDa

Observed MW:

28-35 kDa

Purification Method:

Protein G purification

CloneNo.:

1B1G3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68550-1-Ig (MSX2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68550-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 68550-1-Ig (MSX2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68550-1-PBS in a different storage buffer formulation.