

For Research Use Only

MAP2K3 Monoclonal antibody, PBS Only



Catalog Number: 68574-1-PBS

Featured Product

Basic Information

Catalog Number:

68574-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4861

GenBank Accession Number:

BC032478

GeneID (NCBI):

5606

UNIPROT ID:

P46734

Full Name:

mitogen-activated protein kinase kinase 3

Calculated MW:

347 aa, 39 kDa

Observed MW:

36-40 kDa

Purification Method:

Protein G purification

CloneNo.:

1H3E10

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68574-1-Ig (MAP2K3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68574-1-PBS in a different storage buffer formulation.



WB result of MAP2K3 antibody (68574-1-Ig; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-MAP2K3 transfected A431 cells. This data was developed using the same antibody clone with 68574-1-PBS in a different storage buffer formulation.