

For Research Use Only

SCNN1A Monoclonal antibody, PBS Only



Catalog Number: 68581-1-PBS

Basic Information

Catalog Number: 68581-1-PBS	GenBank Accession Number: BC006526	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6337	CloneNo.: 1C2C2
Source: Mouse	UNIPROT ID: P37088	
Isotype: IgG2a	Full Name: sodium channel, nonvoltage-gated 1 alpha	
Immunogen Catalog Number: AG29548	Calculated MW: 76 kDa	
	Observed MW: ~80 kDa	

Applications

Tested Applications:
WB, Indirect ELISA
Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68581-1-Ig (SCNN1A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68581-1-PBS in a different storage buffer formulation.