

For Research Use Only

TMEM132A Monoclonal antibody, PBS Only



Catalog Number: 68589-1-PBS

Basic Information

Catalog Number: 68589-1-PBS	GenBank Accession Number: BC114569	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 54972	CloneNo.: 1H7G8
Source: Mouse	UNIPROT ID: Q24JP5	
Isotype: IgG2a	Full Name: transmembrane protein 132A	
Immunogen Catalog Number: AG18149	Calculated MW: 1023 aa, 110 kDa	
	Observed MW: 100-140 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68589-1-Ig (TMEM132A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68589-1-PBS in a different storage buffer formulation.