

For Research Use Only

TCF4 Monoclonal antibody, PBS Only

Catalog Number: 68607-1-PBS



Basic Information

Catalog Number: 68607-1-PBS	GenBank Accession Number: BC125084	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6925	CloneNo.: 1C9A1
Source: Mouse	UNIPROT ID: P15884	
Isotype: IgG1	Full Name: transcription factor 4	
Immunogen Catalog Number: AG17897	Calculated MW: 671 aa, 72 kDa	
	Observed MW: 72 kDa; 68 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, rat

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68607-1-Ig (TCF4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68607-1-PBS in a different storage buffer formulation.