

For Research Use Only

# HGF Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68660-1-PBS



## Basic Information

Catalog Number:

68660-1-PBS

Size:

100ug, Concentration: 1mg/ml by  
Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

HZ-1084

GenBank Accession Number:

GeneID (NCBI):

3082

Full Name:

hepatocyte growth factor  
(hepapoietin A; scatter factor)

Purification Method:

Protein A purification

CloneNo.:

2G2C3

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Product Information

68660-1-PBS targets HGF as part of a matched antibody pair:

MP50034-1: 68660-1-PBS capture and 68660-2-PBS detection (validated in Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP50034-1, HGF Monoclonal Matched Antibody Pair - PBS only. 68660-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard HZ-1084. 68660-2-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL