

For Research Use Only

# FLT3 Ligand Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68664-2-PBS



## Basic Information

Catalog Number:

68664-2-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1151

GenBank Accession Number:

GeneID (NCBI):

2323

Full Name:

fms-related tyrosine kinase 3 ligand

Purification Method:

Protein G purification

CloneNo.:

4C4G8

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Product Information

68664-2-PBS targets FLT3 Ligand as part of a matched antibody pair:

MP50046-1: 68664-1-PBS capture and 68664-2-PBS detection (validated in Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50046-1, FLT3LG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68664-1-PBS. Detection antibody: 68664-2-PBS. Standard: HZ-1151. Range: 15.6-1000 pg/mL.