

For Research Use Only

IL-17F Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68666-2-PBS



Basic Information

Catalog Number:

68666-2-PBS

Size:

100ug, Concentration: 1mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1116

GenBank Accession Number:

GeneID (NCBI):

112744

Full Name:

interleukin 17F

Purification Method:

Protein G purification

CloneNo.:

1C5C6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Product Information

68666-2-PBS targets IL-17F as part of a matched antibody pair:

MP50006-1: 68666-1-PBS capture and 68666-2-PBS detection (validated in Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50006-1, IL-17F Monoclonal Matched Antibody Pair, PBS Only. 68666-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard HZ-1116. 68666-2-PBS was HRP conjugated as the detection antibody Range: 62.5-4000 pg/mL