

For Research Use Only

IL-1RL1/ST2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68680-1-PBS



Basic Information

Catalog Number: 68680-1-PBS	GenBank Accession Number: BC030975	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9173	CloneNo.: 4G4B10
Source: Mouse	UNIPROT ID: Q01638-2	
Isotype: IgG1	Full Name: interleukin 1 receptor-like 1	
	Calculated MW: 63 kDa, 37 kDa, 30 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

68680-1-PBS targets IL-1RL1/ST2 as part of a matched antibody pair:

MP50068-1: 68680-1-PBS capture and 68680-2-PBS detection (validated in Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50068-1, IL1RL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68680-1-PBS. Detection antibody: 68680-2-PBS. Standard: Eg0116. Range: 31.25-2000 pg/mL