

For Research Use Only

TNFSF18 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68684-1-PBS



Basic Information

Catalog Number: 68684-1-PBS	GenBank Accession Number: BC069319	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 8995	CloneNo.: 1G11F11
Source: Mouse	UNIPROT ID: Q9UNG2	
Isotype: IgG1	Full Name: tumor necrosis factor (ligand) superfamily, member 18	
	Calculated MW: 177 aa, 20 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Product Information

68684-1-PBS targets TNFSF18 as part of a matched antibody pair:

MP50040-1: 68684-1-PBS capture and 68684-2-PBS detection (validated in Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50040-1, TNFSF18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68684-1-PBS. Detection antibody: 68684-2-PBS. Standard: Eg32094. Range: 62.5-4000 pg/mL.