

For Research Use Only

PDGF-D Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68701-2-PBS



Basic Information

Catalog Number:

68701-2-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29199

GenBank Accession Number:

BC030645

GeneID (NCBI):

80310

UNIPROT ID:

Q9GZPO

Full Name:

platelet derived growth factor D

Calculated MW:

43 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G9C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

68701-2-PBS targets PDGF-D as part of a matched antibody pair:

MP50964-1: 68701-2-PBS capture and 68701-1-PBS detection (validated in Cytometric bead array)

MP50964-4: 68701-2-PBS capture and 68701-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50964-1, PDGFD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68701-2-PBS. Detection antibody: 68701-1-PBS. Standard: Ag5154. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50964-4, PDGFD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68701-2-PBS. Detection antibody: 68701-5-PBS. Standard: Ag5154. Range: 0.098-100 ng/mL