

For Research Use Only

THOC3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68812-1-PBS

Basic Information

Catalog Number: 68812-1-PBS	GenBank Accession Number: BC006849	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 84321	CloneNo.: 3D12D5
Source: Mouse	UNIPROT ID: Q96J01	
Isotype: IgG1	Full Name: THO complex 3	
Immunogen Catalog Number: AG35578	Calculated MW: 327 aa, 36 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

68812-1-PBS targets THOC3 as part of a matched antibody pair:

MP50185-1: 68812-1-PBS capture and 68812-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50185-1, THOC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68812-1-PBS. Detection antibody: 68812-2-PBS. Standard: Ag35578. Range: 0.391-100 ng/mL.