

For Research Use Only

ZC3H12D Monoclonal antibody, PBS Only



Catalog Number: 68847-1-PBS

Basic Information

Catalog Number:

68847-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG21770

GenBank Accession Number:

BC157832

GeneID (NCBI):

340152

UNIPROT ID:

A2A288

Full Name:

zinc finger CCCH-type containing 12D

Calculated MW:

527 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

1E2A2

Applications

Tested Applications:

ELISA, WB

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68847-1-Ig (ZC3H12D antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68847-1-PBS in a different storage buffer formulation.