

For Research Use Only

GPC1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68849-1-PBS

Basic Information

Catalog Number: 68849-1-PBS	GenBank Accession Number: BC051279	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2817	CloneNo.: 2D9E7
Source: Mouse	UNIPROT ID: P35052	
Isotype: IgG1	Full Name: glypican 1	
	Calculated MW: 558 aa, 62 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

68849-1-PBS targets GPC1 as part of a matched antibody pair:

MP50239-1: 68849-1-PBS capture and 68849-2-PBS detection (validated in Cytometric bead array)

MP50239-2: 68849-1-PBS capture and 68849-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50239-1, GPC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68849-1-PBS. Detection antibody: 68849-2-PBS. Standard:Eg0737. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50239-2, GPC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68849-1-PBS. Detection antibody: 68849-3-PBS. Standard:Eg0737. Range: 0.781-100 ng/mL.