

For Research Use Only

CHMP1A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68862-2-PBS

Basic Information

Catalog Number: 68862-2-PBS	GenBank Accession Number: BC007527	Purification Method: Protein A Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5119	CloneNo.: 1D7E7
Source: Mouse	UNIPROT ID: Q9HD42	
Isotype: IgG2a	Full Name: chromatin modifying protein 1A	
Immunogen Catalog Number: AG8376	Calculated MW: 196 aa, 22 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

68862-2-PBS targets CHMP1A as part of a matched antibody pair:

MP50256-1: 68862-1-PBS capture and 68862-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50256-1, CHMP1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68862-1-PBS. Detection antibody: 68862-2-PBS. Standard: Ag8376. Range: 0.781-100 ng/mL.