

For Research Use Only

MDH2 Monoclonal antibody, PBS Only



Catalog Number: 68878-1-PBS

Basic Information

Catalog Number: 68878-1-PBS	GenBank Accession Number: BC001917	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4191	CloneNo.: 2A9G11
Source: Mouse	UNIPROT ID: P40926	
Isotype: IgG2b	Full Name: malate dehydrogenase 2, NAD (mitochondrial)	
Immunogen Catalog Number: AG7750	Calculated MW: 36 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68878-1-Ig (MDH2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68878-1-PBS in a different storage buffer formulation.