

For Research Use Only

CENPJ Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68881-2-PBS

Basic Information

Catalog Number: 68881-2-PBS	GenBank Accession Number: BC024209	Purification Method: Protein A Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 55835	CloneNo.: 1E9B9
Source: Mouse	UNIPROT ID: Q9HC77	
Isotype: IgG2a	Full Name: centromere protein J	
Immunogen Catalog Number: AG28390	Calculated MW: 1338 aa, 153 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

68881-2-PBS targets CENPJ as part of a matched antibody pair.

MP50279-1: 68881-1-PBS capture and 68881-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50279-1, CENPJ Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68881-1-PBS. Detection antibody: 68881-2-PBS. Standard: Ag28390. Range: 6.25-200 ng/mL