

For Research Use Only

# PAX2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68885-1-PBS

## Basic Information

<b>Catalog Number:</b> 68885-1-PBS	<b>GenBank Accession Number:</b> BC148710	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 100ug, Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 5076	<b>CloneNo.:</b> 2E6G6
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q02962	
<b>Isotype:</b> IgG1	<b>Full Name:</b> paired box 2	
<b>Immunogen Catalog Number:</b> AG28286	<b>Calculated MW:</b> 431 aa, 46 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

68885-1-PBS targets PAX2 as part of a matched antibody pair:

MP50289-1: 68885-1-PBS capture and 68886-1-PBS detection (validated in Cytometric bead array)

MP50289-3: 68885-1-PBS capture and 68885-4-PBS detection (validated in Cytometric bead array)

MP50289-4: 68885-1-PBS capture and 68885-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Cytometric bead array standard curve of MP50289-1, PAX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68885-1-PBS. Detection antibody: 68886-1-PBS. Standard:Ag28286. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50289-3, PAX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68885-1-PBS. Detection antibody: 68885-4-PBS. Standard:Ag28286. Range: 6.25-100 ng/mL.



Cytometric bead array standard curve of MP50289-4, PAX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68885-1-PBS. Detection antibody: 68885-5-PBS. Standard:Ag28286. Range: 1.563-100 ng/mL.