

For Research Use Only

# WT1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68912-1-PBS

## Basic Information

<b>Catalog Number:</b> 68912-1-PBS	<b>GenBank Accession Number:</b> BC032861	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 7490	<b>CloneNo.:</b> 1B10D6
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P19544	
<b>Isotype:</b> IgG1	<b>Full Name:</b> Wilms tumor 1	
<b>Immunogen Catalog Number:</b> AG19541	<b>Calculated MW:</b> 449 aa, 49 kDa, 57 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

68912-1-PBS targets WT1 as part of a matched antibody pair.

MP50330-1: 68912-1-PBS capture and 68912-2-PBS detection (validated in Cytometric bead array)

MP50330-2: 68912-1-PBS capture and 68912-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50330-1, WT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68912-1-PBS. Detection antibody: 68912-2-PBS. Standard: Ag19541. Range: 6.25-200 ng/mL



Cytometric bead array standard curve of MP50330-2, WT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68912-1-PBS. Detection antibody: 68912-3-PBS. Standard: Ag19541. Range: 6.25-200 ng/mL