

For Research Use Only

# PARK7/DJ-1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68915-5-PBS

## Basic Information

**Catalog Number:**

68915-5-PBS

**Size:**

100ug, Concentration: 1 mg/ml by  
Nanodrop;

**Source:**

Mouse

**Isotype:**

IgG1

**Immunogen Catalog Number:**

AG28526

**GenBank Accession Number:**

BC008188

**GeneID (NCBI):**

11315

**UNIPROT ID:**

Q99497

**Full Name:**

Parkinson disease (autosomal  
recessive, early onset) 7

**Calculated MW:**

189 aa, 20 kDa

**Purification Method:**

Protein G Magarose purification

**CloneNo.:**

4G9H9

## Applications

**Tested Applications:**

Cytometric bead array, Indirect ELISA

**Species Specificity:**

human

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50337-3, PARK7/DJ-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68915-4-PBS. Detection antibody: 68915-5-PBS. Standard: Ag28526. Range: 1.563-100 ng/mL.