

For Research Use Only

DNMT3B Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68919-6-PBS

Basic Information

Catalog Number: 68919-6-PBS	GenBank Accession Number: NM_001207055	Purification Method: Protein A Magarose purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1789	CloneNo.: 1F12B3
Source: Mouse	UNIPROT ID: Q9UBC3	
Isotype: IgG2b	Full Name: DNA (cytosine-5-)-methyltransferase 3 beta	
Immunogen Catalog Number: AG33907	Calculated MW: 95 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data