

For Research Use Only

GPR64 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68963-2-PBS



Basic Information

Catalog Number:

68963-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31274

GenBank Accession Number:

BC099901

GeneID (NCBI):

10149

UNIPROT ID:

Q8LZP9

Full Name:

G protein-coupled receptor 64

Calculated MW:

1017 aa, 112 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D10C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

68963-2-PBS targets GPR64 as part of a matched antibody pair.

MP50415-1: 68963-1-PBS capture and 68963-2-PBS detection (validated in Cytometric bead array)

MP50415-2: 68963-3-PBS capture and 68963-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50415-1, GPR64 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68963-1-PBS. Detection antibody: 68963-2-PBS. Standard: Ag31274. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50415-2, GPR64 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68963-3-PBS. Detection antibody: 68963-2-PBS. Standard: Ag31274. Range: 0.195-100 ng/mL