

For Research Use Only

CMPK2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68974-3-PBS



Basic Information

Catalog Number: 68974-3-PBS	GenBank Accession Number: BC141802	Purification Method: Protein A Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 129607	CloneNo.: 2C9G12
Source: Mouse	UNIPROT ID: Q5EBMO	
Isotype: IgG2b	Full Name: cytidine monophosphate (UMP-CMP) kinase 2, mitochondrial	
Immunogen Catalog Number: AG22885	Calculated MW: 449 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

68974-3-PBS targets CMPK2 as part of a matched antibody pair.

MP50436-2: 68974-1-PBS capture and 68974-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50436-2, CMFK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68974-1-PBS. Detection antibody: 68974-3-PBS. Standard: Ag22885. Range: 0.391-50 ng/mL