

For Research Use Only

IL-10 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 69018-2-PBS



Basic Information

Catalog Number: 69018-2-PBS	GenBank Accession Number: GeneID (NCBI): 3586	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	Full Name: interleukin 10	CloneNo.: 7B8G3
Source: Mouse		
Isotype: IgG2b		
Immunogen Catalog Number: HZ-1145		

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

69018-2-PBS targets IL-10 as part of a matched antibody pair.

MP50035-1: 69018-1-PBS capture and 69018-2-PBS detection (validated in Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50035-1, IL-10
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 69018-1-PBS. Detection
antibody: 69018-2-PBS. Standard: HZ-1145. Range:
15.6-1000 pg/mL.