

For Research Use Only

AMPK alpha Recombinant antibody, PBS Only



Catalog Number: 80209-5-PBS

Basic Information

Catalog Number:

80209-5-PBS

Size:

100ug, Concentration: 1mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1367

GenBank Accession Number:

BC012622

GeneID (NCBI):

5562

UNIPROT ID:

Q13131

Full Name:

protein kinase, AMP-activated, alpha
1 catalytic subunit

Calculated MW:

63 kDa

Purification Method:

Protein A purification

CloneNo.:

4N5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 80209-5-RR (AMPK alpha antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 80209-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 80209-5-RR against Human AMPK alpha were performed. The affinity constant is below 1 pM.