

For Research Use Only

SNAI1 Recombinant antibody, PBS Only



Catalog Number: 81584-5-PBS

Basic Information

Catalog Number: 81584-5-PBS	GenBank Accession Number: BC012910	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6615	CloneNo.: 240164C11
Source: Rabbit	UNIPROT ID: O95863	
Isotype: IgG	Full Name: snail homolog 1 (Drosophila)	
Immunogen Catalog Number: AG24248	Calculated MW: 264 aa, 29 kDa	

Applications

Tested Applications:
ELISA, WB

Species Specificity:
human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 81584-5-RR (SNAI1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 81584-5-PBS in a different storage buffer formulation.



1x10⁶ A549 cells were intracellularly stained with 0.25 ug SNAI1 Recombinant antibody (81584-5-RR, Clone:240164C11) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with True-Nuclear Transcription Factor Buffer Set. This data was developed using the same antibody clone with 81584-5-PBS in a different storage buffer formulation.