

For Research Use Only

IL1B Recombinant antibody

Catalog Number: 82696-5-RR



Basic Information

| | | |
|--------------------------------------|--|---|
| Catalog Number: 82696-5-RR | GenBank Accession Number: BC008678 | Purification Method: Protein A purification |
| Size: 100ul | GeneID (NCBI): 3553 | CloneNo.: 2E22 |
| Source: Rabbit | ENSEMBL Gene ID: ENSG00000125538 | |
| Isotype: IgG | UNIPROT ID: P01584 | |
| | Full Name: interleukin 1, beta | |
| | Calculated MW: 269 aa, 31 kDa | |

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data