

For Research Use Only

IL1B Recombinant antibody

Catalog Number: 82696-8-RR



Basic Information

Catalog Number:
82696-8-RR

Size:
100ul

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
BC008678

GeneID (NCBI):
3553

ENSEMBL Gene ID:
ENSG00000125538

UNIPROT ID:
P01584

Full Name:
interleukin 1, beta

Calculated MW:
269 aa, 31 kDa

Purification Method:
Protein A purification

CloneNo.:
2G9

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data