

For Research Use Only

# GRP94 Recombinant antibody, PBS Only



Catalog Number: 82805-1-PBS

## Basic Information

|  |  |   |
|--|--|---|
| <b>Catalog Number:</b><br>82805-1-PBS                      | <b>GenBank Accession Number:</b><br>BC066656                         | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>100ug , Concentration: 1mg/ml by Nanodrop; | <b>GeneID (NCBI):</b><br>7184  | <b>CloneNo.:</b><br>6K5                               |
| <b>Source:</b><br>Rabbit                                   | <b>UNIPROT ID:</b><br>P14625   |   |
| <b>Isotype:</b><br>IgG                                     | <b>Full Name:</b><br>heat shock protein 90kDa beta (Grp94), member 1 |   |
| <b>Immunogen Catalog Number:</b><br>AG6344                 | <b>Calculated MW:</b><br>92 kDa                                      |   |
|  | <b>Observed MW:</b><br>94-100 kDa                                    |   |

## Applications

**Tested Applications:**  
WB, Indirect ELISA

**Species Specificity:**  
Human

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82805-1-RR (GRP94 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82805-1-PBS in a different storage buffer formulation.