

For Research Use Only

Mouse Snf1lk Recombinant antibody, PBS Only

Catalog Number:82851-3-PBS



Basic Information

Catalog Number:

82851-3-PBS

Size:

100ug , Concentration: 1mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33435

GenBank Accession Number:

NM_010831

GeneID (NCBI):

17691

UNIPROT ID:

Q60670

Full Name:

salt inducible kinase 1

Calculated MW:

85 kDa

Observed MW:

85 kDa

Purification Method:

Protein A purification

CloneNo.:

3D2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82851-3-RR (Snf1lk antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82851-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82851-3-RR against Mouse Snf1lk were performed. The affinity constant is 0.686 nM.