

For Research Use Only

SFTPB Recombinant antibody, PBS Only



Catalog Number: 82866-16-PBS

Basic Information

Catalog Number:

82866-16-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23605

GenBank Accession Number:

BC032785

GeneID (NCBI):

6439

UNIPROT ID:

P07988

Full Name:

surfactant protein B

Calculated MW:

381 aa, 42 kDa

Purification Method:

Protein A purification

CloneNo.:

230240D4

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 82866-16-RR (SFTPb antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82866-16-PBS in a different storage buffer formulation.