

For Research Use Only

TBCD Recombinant antibody, PBS Only

Catalog Number: 82911-2-PBS



Basic Information

Catalog Number: 82911-2-PBS	GenBank Accession Number: BC003094	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6904	CloneNo.: 230209D7
Source: Rabbit	UNIPROT ID: Q9BTW9	
Isotype: IgG	Full Name: tubulin folding cofactor D	
Immunogen Catalog Number: AG6310	Calculated MW: 133 kDa	

Applications

Tested Applications:
Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 82911-2-RR against Human TBCD were performed. The affinity constant is 2.07 nM.