

For Research Use Only

MYOF Recombinant antibody, PBS Only



Catalog Number: 82917-7-PBS

Basic Information

Catalog Number:

82917-7-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_133337

GeneID (NCBI):

26509

UNIPROT ID:

Q9NZM1

Full Name:

myoferlin

Calculated MW:

235 kDa

Purification Method:

Protein A purification

CloneNo.:

230213B6

Applications

Tested Applications:

IF, ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using MYOF antibody (82917-7-RR, Clone: 230213B6) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1). This data was developed using the same antibody clone with 82917-7-PBS in a different storage buffer formulation.