

For Research Use Only

NFX1 Recombinant antibody, PBS Only



Catalog Number: 82924-2-PBS

Basic Information

Catalog Number:

82924-2-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9448

GenBank Accession Number:

BC012151

GeneID (NCBI):

4799

UNIPROT ID:

Q12986

Full Name:

nuclear transcription factor, X-box binding 1

Calculated MW:

1120 aa, 124 kDa

Observed MW:

124 kDa

Purification Method:

Protein A purification

CloneNo.:

230211H4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82924-2-RR (NFX1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82924-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82924-2-RR against Human NFX1 were performed. The affinity constant is below 1 pM.