

For Research Use Only

# IFT140 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82947-6-PBS

## Basic Information

<b>Catalog Number:</b> 82947-6-PBS	<b>GenBank Accession Number:</b> BC035577	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 9742	<b>CloneNo.:</b> 230231H2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q96RY7	
<b>Isotype:</b> IgG	<b>Full Name:</b> intraflagellar transport 140 homolog (Chlamydomonas)	
<b>Immunogen Catalog Number:</b> AG11119	<b>Calculated MW:</b> 1462 aa, 165 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Product Information

82947-6-PBS targets IFT140 as part of a matched antibody pair:

MP00029-2: 82947-4-PBS capture and 82947-6-PBS detection (validated in Cytometric bead array)

MP00029-3: 82947-3-PBS capture and 82947-6-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
100% PBS pH 7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00029-2, IFT140 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82947-4-PBS. Detection antibody: 82947-6-PBS. Standard: Ag11119. Range: 2.5-20 ng/mL



Cytometric bead array standard curve of MP00029-3, IFT140 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82947-3-PBS. Detection antibody: 82947-6-PBS. Standard: Ag11119. Range: 2.5-20 ng/mL