

For Research Use Only

CPO Recombinant antibody, PBS Only



Catalog Number: 82952-2-PBS

Basic Information

Catalog Number: 82952-2-PBS	GenBank Accession Number: BC112078	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 130749	CloneNo.: 230153F4
Source: Rabbit	UNIPROT ID: Q8IVL8	
Isotype: IgG	Full Name: carboxypeptidase O	
Immunogen Catalog Number: AG24059	Calculated MW: 374 aa, 43 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 82952-2-RR (CPO antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82952-2-PBS in a different storage buffer formulation.