

For Research Use Only

AFP Recombinant antibody, PBS Only (Detector)

Catalog Number: 82982-2-PBS



Basic Information

Catalog Number: 82982-2-PBS	GenBank Accession Number: BC027881	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 174	CloneNo.: 230435F11
Source: Rabbit	UNIPROT ID: P02771	
Isotype: IgG	Full Name: alpha-fetoprotein	
	Calculated MW: 69 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

82982-2-PBS targets AFP as part of a matched antibody pair.

MP00001-1: 82982-1-PBS capture and 82982-2-PBS detection (validated in Cytometric bead array)

MP00001-2: 82982-3-PBS capture and 82982-2-PBS detection (validated in Cytometric bead array)

MP00001-3: 82982-4-PBS capture and 82982-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00001-1, AFP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82982-1-PBS. Detection antibody: 82982-2-PBS. Standard: Eg0539. Range: 1.25-20 ng/mL



Cytometric bead array standard curve of MP00001-2, AFP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82982-3-PBS. Detection antibody: 82982-2-PBS. Standard: Eg0539. Range: 1.25-20 ng/mL



Cytometric bead array standard curve of MP00001-3, AFP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82982-4-PBS. Detection antibody: 82982-2-PBS. Standard: Eg0539. Range: 0.625-20 ng/mL