

For Research Use Only

ANO9 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83025-3-PBS

Basic Information

Catalog Number: 83025-3-PBS	GenBank Accession Number: NM_001012302	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 338440	CloneNo.: 230256D2
Source: Rabbit	UNIPROT ID: A1A5B4	
Isotype: IgG	Full Name: anoctamin 9	
	Calculated MW: 90 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83025-3-PBS targets ANO9 as part of a matched antibody pair.

MP00092-2: 83025-3-PBS capture and 83025-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00092-2, ANO9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83025-3-PBS. Detection antibody: 83025-4-PBS. Standard: SY00332. Range: 0.078-5 ng/mL