

For Research Use Only

FGF21 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83033-1-PBS



Basic Information

Catalog Number: 83033-1-PBS	GenBank Accession Number: BC018404	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 26291	CloneNo.: 230236B2
Source: Rabbit	UNIPROT ID: Q9NSA1	Affinity: $K_D = 3.67 \times 10^{-10}M$
Isotype: IgG	Full Name: fibroblast growth factor 21	$K_{Off} = 3.38 \times 10^{-4}M$
	Calculated MW: 22 kDa	$K_{On} = 9.20 \times 10^5M$

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Product Information

83033-1-PBS targets FGF21 as part of a matched antibody pair:

MP00005-1: 83033-1-PBS capture and 83033-2-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00005-1, FGF21 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83033-1-PBS. Detection antibody: 83033-2-PBS. Standard: Eg0102. Range: 31.25-2000 pg/mL